



Only **ACUVUE®** has
BLINK STABILIZED®
Design

for exceptional vision,
stability, and comfort¹



SCAN for
FITTING CALCULATOR

up to
-2.75D

around the clock
coverage²

89%

coverage for
astigmatic patients²

99%

first lens fit
success³

Current SKUs
 Expanded Coverage SKUs
 Available in 30-packs only

ACUVUE® OASYS
1-Day for
ASTIGMATISM
Contact Lens
Parameters

		LOW MINUS					HIGH MINUS					LOW PLUS					HIGH PLUS				
		Plano to -6.00D in 0.25D steps					-6.50D to -9.00D in 0.50D steps					+0.25D to +4.00D in 0.25D steps					+4.25D to +6.00D in 0.25D steps				
		-0.75	-1.25	-1.75	-2.25	-2.75	-0.75	-1.25	-1.75	-2.25	-2.75	-0.75	-1.25	-1.75	-2.25	-2.75	-0.75	-1.25	-1.75	-2.25	-2.75
WITH THE RULE	10°																				
	20°																				
	160°																				
	170°																				
	180°																				
AGAINST THE RULE	70°																				
	80°																				
	90°																				
	100°																				
	110°																				
OBLIQUES	30°																				
	40°																				
	50°																				
	60°																				
	120°																				
	130°																				
	140°																				
	150°																				

^Rx range: -1.50 to -4.00DS inclusive, -0.75 to 1.50DC inclusive and axes within 180 ± 15 degrees or 90 ± 15 degrees.

References: 1. JJV Data on file 2024, ACUVUE® Brand Contact Lenses for ASTIGMATISM overall fitting success, orientation position, rotational stability and vision performance. 2. JJV Data on file. SKU Coverage Data.
 3. Straker, B, Hamada, W, Sulley, A, et al. Fitting performance and efficiency with a new silicone hydrogel daily disposable toric contact lens. Contact Lens Anterior Eye. 2018;41(S1):S57-S58. doi: 10.1016/j.clae.2018.03.049.

Important Safety Information: ACUVUE® Contact Lenses are indicated for vision correction. As with any contact lens, eye problems, including corneal ulcers, can develop. Some wearers may experience mild irritation, itching or discomfort. Lenses should not be prescribed if patients have any eye infection, or experience eye discomfort, excessive tearing, vision changes, redness or other eye problems. Consult the package insert for complete information. Complete information is also available from Johnson & Johnson Vision by calling 1-800-267-5098 or by visiting www.jnjvisionpro.com/en-ca/.

© Johnson & Johnson and its affiliates 2026 2026PP11182